

50HZ

Ambient	Model	Heat Generated W 30°C→-86°C	Heat Generated W/Hr @-86°C	Power Consumption W 30°C→-86°C	Power Consumption W/Hr @-86°C	Power Consumption W/Month @-86°C
32°C	AV039	7250	608	3470	490	14700
	AV050	8700	730	4170	520	15600
	AV061	9152	826	4490	690	20700
	AV073	10023	910	4540	760	22800
25°C	AV039	6524	550	3120	440	13200
	AV050	7920	660	3750	470	14100
	AV061	8242	669	4040	620	18600
	AV073	9060	745	4430	690	20700